



PRICE LIST



TEST CODE	TEST NAME	TEST PRICE	TEST SPECIMEN
HEMATOLOGY & COAGULATION			
PV001	HEMOGLOBIN; HB	150	3 mL (1 mL min) Whole Blood in 1 Lavender Top (EDTA) Ship refrigerated. DO NOT FREEZE
PV008	COMPLETE BLOOD COUNT; CBC	400	3 mL (1 mL min) Whole Blood in 1 Lavender Top (EDTA) Ship refrigerated. DO NOT FREEZE
PV009	HEMOGRAM *CBC * ESR	450	3 mL (1 mL min) Whole Blood in 1 Lavender Top (EDTA) Ship refrigerated. DO NOT FREEZE
PV011	DLC (DIFFERENTIAL LEUCOCYTE COUNT)	150	3 mL (1 mL min) Whole Blood in 1 Lavender Top (EDTA) Ship refrigerated. DO NOT FREEZE
PV012	TLC (TOTAL LEUCOCYTE COUNT)	170	3 mL (1 mL min) Whole Blood in 1 Lavender Top (EDTA) Ship refrigerated. DO NOT FREEZE
PV018	PLATELET COUNT	150	3 mL (1 mL min) Whole Blood in 1 Lavender Top (EDTA) Ship refrigerated. DO NOT FREEZE
PV111	BLOOD SMEAR -HEAMOPROTOZOAN	300	Peripheral whole blood or 2 Peripheral blood smears. Give brief clinical history.
PV091	KNOTT'S HEARTWORM TEST	400	3 mL (1 mL min) Whole Blood in 1 Lavender Top (EDTA) Ship refrigerated. DO NOT FREEZE
PV010	PROTHROMBIN TIME STUDIES	500	3 mL (2 mL min) Whole Blood in 1 Blue Top (Sodium Citrate) tube with proper dilution
PV025	PARTIAL THROMBOPLASTIN TIME, ACTIVATED; PTT; APTT	600	3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3000 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3000 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred.
PV046	IRON STUDIES (SERUM IRON, TOTAL IRON BINDING CAPACITY & TRANSFERRIN SATURATION	1200	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV211	RETICULOCYTE COUNT	350	3 mL (1 mL min) Whole Blood in 1 Lavender Top (EDTA)
BIOCHEMISTRY			
PV020	PROTEIN, TOTAL, SERUM *TOTAL PROTEIN *ALBUMIN *A : G RATIO	150	2 mL (1 mL min) serum from 1 SST Ship refrigerated or frozen.
PV021	ALBUMIN, SERUM	130	2 mL (1 mL min) serum from 1 SST Ship refrigerated or frozen.
PV024	AMYLASE, SERUM	400	2 mL (1 mL min) serum from 1 SST Ship refrigerated or frozen.
PV042	LIPASE, SERUM	580	2 mL (1 mL min) serum from 1 SST Ship refrigerated or frozen.
PV026	BICARBONATE	550	2 mL (1 mL min) serum from 1 SST Ship refrigerated or frozen.
PV027	BILIRUBIN, DIRECT	140	2 mL (1 mL min) serum from 1 SST. Wrap tube in aluminium foil to protect from light.
PV028	BILIRUBIN TOTAL	150	2 mL (1 mL min) serum from 1 SST. Wrap tube in aluminium foil to protect from light.

TEST CODE	TEST NAME	TEST PRICE	TEST SPECIMEN
BIOCHEMISTRY			
PV110	BILIRUBIN TOTAL, DIRECT & INDIRECT	250	2 mL (1 mL min.) serum from 1 SST. Wrap tube in aluminium foil to protect from light. Ship refrigerated or frozen.
PV022	ALKALINE PHOSPHATASE; ALP	150	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV023	AMMONIA, BLOOD	950	3 mL (1 mL min) Whole Blood in 1 Lavender Top (EDTA) tube separate IMMEDIATELY, Ship frozen, AVOID hemolysis
PV047	SGOT	150	2 mL (1 mL min) serum from 1 SST Ship refrigerated or frozen.
PV048	SGPT	150	2 mL (1 mL min) serum from 1 SST Ship refrigerated or frozen.
PV035	GGTP; GGT; GAMMA GT; GAMMA GLUTAMYL TRANSPEPTIDASE	250	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV036	GLUCOSE FASTING (F)	120	2 mL (1 mL min) plasma from 1 Grey Top (Sodium Fluoride) tube. Ship refrigerated or frozen.
PV037	GLUCOSE, RANDOM	120	2 mL (1 mL min) plasma from 1 Grey Top (Sodium Fluoride) tube. Ship refrigerated or frozen.
PV038	GLUCOSE, POSTPRANDIAL (PP), 2 HOURS	120	2 mL (1 mL min) plasma from 1 Grey Top (Sodium Fluoride) tube. Ship refrigerated or frozen.
PV032	CHOLESTEROL, TOTAL	180	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV040	CHOLESTEROL, HDL, DIRECT	220	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV041	CHOLESTEROL, LDL, DIRECT	220	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV050	TRIGLYCERIDES, SERUM	260	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV043	MAGNESIUM, SERUM	450	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV030	CALCIUM, SERUM	150	2 mL (1 mL min) serum from 1 SST. Avoid use of the tourniquet.
PV031	CHLORIDE, SERUM	150	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV044	PHOSPHORUS, SERUM	150	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV045	POTASSIUM, SERUM	150	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV049	SODIUM, SERUM	150	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV051	UREA , SERUM	150	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV029	BUN; BLOOD UREA NITROGEN	150	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV034	CREATININE, SERUM	150	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV052	URIC ACID, SERUM	150	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.

TEST CODE	TEST NAME	TEST PRICE	TEST SPECIMEN
BIOCHEMISTRY			
PV130	Vitamin D, 25 - Hydroxy	1500	2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.
PV089	VITAMIN B12; CYANOCOBALAMIN	1300	3 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Wrap in aluminium foil to protect from light. Overnight fasting is mandatory.
PV099	LIPID PROFILE, BASIC	880	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Minimum 12 hours fasting will give good results. Grossly hemolysed samples are not acceptable.
PV101	ELECTROLYTES, SERUM *SODIUM * POTASSIUM * CHLORIDE	360	2 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.
PV102	KIDNEY PANEL; KFT	700	3 mL (1.5 mL min.) serum from 1 SST. Overnight fasting is mandatory. Ship refrigerated or frozen.
PV103	LIVER PANEL 1; LFT	700	2 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.
PV112	PHENOBARBITONE	880	2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: just before the next dose is due at steady state
PV146	LDH; LACTIC DEHYDROGENASE	320	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated. DO NOT FREEZE. Hemolysed specimens are not acceptable.
PV033	CPK (CREATININE KINASE)	350	2 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.
PV236	FRUCTOSAMINE	470	2 mL (1 mL min.) serum from 1 SST.
PV242	VITAMIN D 1, 25 DIHYDROXY	3200	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.
CLINICAL PATH & URINE BIOCHEMISTRY			
PV053	URINE EXAMINATION, ROUTINE; URINE, R/E	150	10 mL (5 mL min) aliquot for first morning urine in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE
PV054	SPECIFIC GRAVITY, URINE (By species specific Refractometer)	150	10 mL (5 mL min) aliquot for first morning urine in a sterile screw capped container.
PV055	PROTEIN, TOTAL, RANDOM URINE	150	10 mL (5 mL min) aliquot for first morning urine in a sterile screw capped container.
PV056	URINE, GLUCOSE	150	10 mL (5 mL min) aliquot for first morning urine in a sterile screw capped container.
PV120	MICROALBUMIN 1ST MORNING/RANDOM URINE	500	10 ml (5 ml min.) aliquot of first morning / random urine in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.
PV121	MICROALBUMIN CREATININE RATIO, URINE	500	15 mL (10 mL min.) aliquot of random urine in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.
PV122	PROTEIN CREATININE RATIO, URINE (UP/C)	420	15 mL (10 mL min.) aliquot of first morning / random urine in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.

TEST CODE	TEST NAME	TEST PRICE	TEST SPECIMEN
CLINICAL PATH & URINE BIOCHEMISTRY			
PV080	STOOL EXAMINATION, OCCULT BLOOD	200	5g (3g min) stool in a leak-proof screw capped container. Do not use match box or absorbant material. Outstation specimen not accepted.
PV081	STOOL EXAMINATION, ROUTINE; STOOL, R/E	150	5g (3g min) stool in a leak-proof screw capped container. Do not use match box or absorbant material. Outstation specimen not accepted.
PV082	STOOL EXAMINATION; FLOTATION	200	5g (3g min) stool in a leak-proof screw capped container. Do not use match box or absorbant material. Outstation specimen not accepted.
PV090	SKIN SCRAPPING TEST	400	Best to collect Scrapings in sterile glycerine or liquid paraffin. Sample vial should be tightly packed. Ship refrigerated.
PV075	CYTOLOGY, VAGINAL	350	Prepare the smear immediately after withdrawal of vaginal swab, do not allow smear to air dry before fixation Brief clinical history in Cytology Requisition Form is mandatory.
PV083	SEMEN ANALYSIS; SEMINOGRAM	400	Submit Semen in a sterile screw capped container. Specimen must reach the lab within 1 hour of collection. Print exact time of collection on specimen container and test Requisition form. Outstation samples not accepted
KIT BASED TEST			
PV092	CANINE PARVO VIRUS	700	Faeces Swab/ Serum
PV093	CANINE HEART WORM	700	3 ml Whole Blood, Serum or Plasma
PV094	CANINE DISTEMPER VIRUS	700	Eye mucus conjunctiva or secretions from nasal cavity swab 3 ml Urine, Serum or Plasma
PV095	CANINE LEISHMANIA AB TEST	700	3 ml Whole Blood, Serum or Plasma Ship refrigerated
PV096	CANINE LEPTOSPIROSIS	1200	3 ml Whole Blood, Serum or Plasma Ship refrigerated
PV097	CANINE EHRLICHIA TEST	700	3 ml Whole Blood, Serum or Plasma Ship refrigerated
PV098	FELINE IMMUNODEFICIENCY VIRUS AB, FELINE LEUKEMIA VIRUS AG AND FELINE HEARTWORM AB TEST	2000	3 ml Whole Blood, Serum or Plasma Ship refrigerated
PV209	SNAP 4DX PLUS TEST	2000	3 ml Whole Blood, Serum or Plasma Ship refrigerated
PV219	CANINE PANCREATIC LIPASE (CPL)	2500	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen
PV239	CANINE BRUCELLOSIS AB TEST	800	3 mL (1 mL min.) whole blood EDTA, serum or plasma.
SEROLOGY			
PV123	BRUCELLA ANTIBODY, IgM	1000	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen
PV124	BRUCELLA ANTIBODY, IgG	1000	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.
PV132	IMMUNOGLOBULIN IgM,SERUM	400	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.
PV133	IMMUNOGLOBULIN IgA,SERUM	350	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.

TEST CODE	TEST NAME	TEST PRICE	TEST SPECIMEN
SEROLOGY			
PV152	IMMUNOGLOBULIN IgE, SERUM	690	2 ml (1 ml min.) serum from 1 SST. Ship refrigerated or frozen.
PV134	IMMUNOGLOBULIN IgG, SERUM	350	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.
PV201	BRUCELLA AGGLUTINATION TEST	700	2 ml (1 ml min.) serum from 1 SST Ship refrigerated or frozen.
PV203	LEISHMANIA ANTIBODY, IgG	1600	2 ml (1 ml min.) serum from 1 SST Ship refrigerated or frozen.
CYTOLOGY & HISTOPATHOLOGY			
PV074	CYTOLOGY, FINE NEEDLE ASPIRATION; FNAC	750	Lumps/CT or USG guided deep lumps. Smear on 4 glass slides and fix immediately in absolute alcohol for 30 minutes. Air dry completely. Do not allow smear to air dry before fixation Brief clinical history in Cytology Requisition Form is mandatory.
PV076	CYTOLOGY, BODY FLUIDS FOR MALIGNANT CELLS	500	Slide prepared or body fluids in screw capped container.
PV077	HISTOPATHOLOGY, BIOPSY, SOFT TISSUE (SMALL SPECIMEN)	750	Submit specimen in 10% buffered formalin. Formalin should cover the specimen completely. Ship at room temperature. Brief clinical history in Histology Requisition Form is mandatory.
PV210	HISTOPATHOLOGY, BIOPSY, LARGE SPECIMEN	1190	Submit specimen in 10% buffered formalin. Formalin should cover the specimen completely. Ship at room temperature. Brief clinical history in Histology Requisition Form is mandatory.
PV078	HISTOPATHOLOGY, BIOPSY, BONE	1100	Submit specimen in 10% buffered formalin. Formalin should cover the specimen completely. Ship at room temperature. Brief clinical history in Histology Requisition Form is mandatory.
PV079	HISTOPATHOLOGY, BIOPSY- SECOND OPINION	600	Histopathology section or paraffin embedded tissue blocks previously reported upon.
PV136	HISTOPATHOLOGY, SPECIAL STAIN, MUCICARMINEPERIODIC ACID SCHIFF (PAS)	450	Submit specimen in 10% buffered formalin. Formalin should cover the specimen completely. Ship at room temperature. Brief clinical history in Histology Requisition Form is mandatory.
PV137	HISTOPATHOLOGY, SPECIAL STAIN, MASSON TRICHROME	450	Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide clinical history.
PV073	FLUID EXAMINATION, ROUTINE, BODY FLUIDS INCLUDES ASCITIC/ PERITONEAL/ PERICARDIAL/ PLEURAL/ BRONCHIAL WASHINGS/ PROSTATIC/ SYNOVIAL/ CAPD FLUID DOES NOT INCLUDE PAPANICOLAOU STAIN.	750	5 mL (2 mL min) fluid in a sterile screw capped container. Specify type of fluid on test request form and container.

TEST CODE	TEST NAME	TEST PRICE	TEST SPECIMEN
CYTOLOGY & HISTOPATHOLOGY			
PV228	IMMUNOHISTOCHEMISTRY (IHC)	1570	Submit specimen in 10% buffered formalin. Formalin should cover the specimen completely. Ship at room temperature. Brief clinical history in Histology Requisition Form is mandatory.
PV241	BODY FLUID EXAMINATION ROUTINE (Biochemical, Cytological & Microbiological analysis)	1500	5 ml (3 mL min.) fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE. Specify type of fluid on Test Request Form and container.
MICROBIOLOGY			
PV059	CULTURE, AEROBIC, BLOOD, RAPID INCLUDES ANTIBIOTIC SENSITIVITY	950	Collect 5 mL blood aseptically in a Aerobic blood culture bottle.
PV226	CULTURE, ANAEROBIC, BLOOD, RAPID INCLUDES ANTIBIOTIC SENSITIVITY	1200	Collect 5 mL blood aseptically in a Anaerobic blood culture bottle.
PV060	CULTURE, AEROBIC, BODY FLUIDS, RAPID INCLUDES ANTIBIOTIC SENSITIVITY	1000	Submit aseptically collected body fluids in a sterile vial. Specify collection site and provide brief clinical history
PV064	CULTURE, AEROBIC, MISCELLANEOUS INCLUDES ANTIBIOTIC SENSITIVITY	750	Submit aseptically collected specimen. Specify type of specimen and collection site on test request form. Ship refrigerated.
PV227	CULTURE, ANAEROBIC, MISCELLANEOUS INCLUDES ANTIBIOTIC SENSITIVITY	1200	Submit aseptically collected specimen in anaerobic medium (RCM). Specify type of specimen and collection site on test request form. Ship refrigerated.
PV118	CULTURE, STOOL INCLUDES ANTIBIOTIC SENSITIVITY IF INDICATED.	750	Aseptically collected stool in sterile vial. Ship refrigerated.
PV119	CULTURE, URINE INCLUDES ANTIBIOTIC SENSITIVITY IF INDICATED.	750	Aseptically collected urine in sterile vial. Ship refrigerated.
PV125	BOVINE MILK CULTURE AND SENSITIVITY (Conventional culture and Automated Sensitivity)	500	Aseptically collect milk in sterile vial. Ship refrigerated. Please use separate PV code for individual teat/udder
PV129	MILK CULTURE AND SENSITIVITY (Rapid Automated Culture and Sensitivity)	1000	Aseptically collected milk in sterile vial. Ship refrigerated.
PV071	CULTURED FUNGUS IDENTIFICATION DOES NOT INCLUDE ANTIFUNGAL SENSITIVITY	750	Submit pure isolation appropriate media Submit aseptically collected specimen in a sterile vial.
PV072	CULTURE, FUNGUS INCLUDES ANTIFUNGAL SENSITIVITY FOR YEAST ONLY IF INDICATED	1000	Specify collection site and provide brief clinical history
PV207	FUNGUS EXAMINATION, ROUTINE, KOH CALCOFLUOR WHITE PREPARATION	300	Submit Skin, Nail scrapings, Hair (including root), Cerebrospinal fluid, in a sterile screw capped container OR 2 mL whole blood / bone marrow in 1 Green Top (Sodium Heparin) tube. Ship at 18–22°C. Submit Body fluid, Sputum, Bronchial washings / Tissue in a sterile screw capped container OR Oral/Genital/wound swabs in 0.5 mL sterile normal saline. Ship refrigerated.
HORMONE ANALYSIS			
PV087	T4, TOTAL (CANINE SPECIFIC)	400	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.

TEST CODE	TEST NAME	TEST PRICE	TEST SPECIMEN
HORMONE ANALYSIS			
PV086	T4, FREE; FT4	500	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen
PV115	ACTH	1350	2 mL (0.5 mL min.) plasma from 1 Lavender Top (EDTA) tube. Draw blood between 8-10 a.m. Mix well, separate and transfer plasma IMMEDIATELY into a sterile screw capped vial and FREEZE. Ship frozen. DO NOT THAW.
PV116	CORTISOL	500	2 mL (1 mL min.) serum from 1 sst. Diurnal variation present; take specimen between 7-9 a.m. Specify time and date on sample container and test request form. Ship refrigerated or frozen. Overnight fasting is preferred.
PV154	ESTRADIOL (E2) SERUM	540	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen
PV155	FSH; FOLLICULAR STIMULATING HORMONE	500	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen
PV156	PROLACTIN SERUM	500	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen
PV151	TESTOSTERONE, TOTAL	570	2 mL (1 mL min.) serum from 1 sst. Ship refrigerated or frozen. specify age and sex on test request form.
PV084	PROGESTERONE, SERUM	600	2 mL (1 mL min.) serum from Red Top (avoid serum separator tubes) ship refrigerated or frozen.
PV150	PTH (PARATHYROID HORMONE), INTACT	1500	3 mL (1.5 mL min.) serum from 1 sst drawn between 8-10 a.m in a fasting condition. Allow blood to clot. Separate serum and freeze within 30 minutes. Ship frozen. Do not thaw. Hemolysed samples not acceptable.
SPECIALISED TEST			
PV148	TROPONIN - I	1500	2 mL (1.5 mL min.) serum from 1 SST. ship refrigerated or frozen
PV159	TROPONIN - T	1500	2 mL (1.5 mL min.) serum from 1 SST. ship refrigerated or frozen
PV205	FAVN TEST (RABIES ANTIBODY TITER FOR EXPORT ANIMALS)	16500	5 mL blood in SST or 3 mL Serum collected aseptically
PV204	BABESIA DETECTION, (PCR ARRAY)	1500	3 mL (2 mL min) Whole Blood in 1 Lavender Top (EDTA) Ship refrigerated
PV220	CANINE TICK BORNE DISEASE PCR PANEL	2000	3 mL (2 mL min) Whole Blood in 1 Lavender Top (EDTA) Ship refrigerated
PV221	CANINE EHRLICHIA DETECTION PCR	900	3 mL (2 mL min) Whole Blood in 1 Lavender Top (EDTA) Ship refrigerated
PV145	C-REACTIVE PROTEIN; CRP	350	2 mL (1.5 mL min.) Serum from 1sst. ship refrigerated or frozen Over night fasting is preferred.
PV057	STONE ANALYSIS WITH PICTURE	1500	Submit stone (s) in a clean screw capped container. Ship at room temp. or refrigerated


TEST CODE	TEST NAME	TEST PRICE	TEST SPECIMEN
PANELS			
PV109	LFT*KFT PANELS	1320	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.
PV104	PREANESTHETIC PANEL(CBC*ALP*ALT*TP *ALB*GLU*BUN*CREA*GLOB*A/G)	1250	2 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. 3 mL (1 mL min)Whole Blood in 1 Lavender Top (EDTA) and (Sodium Fluoride) Grey Top tube
PV107	ANEMIA PANEL:COMPLETE BLOOD COUNT * SERUM IRON STUDIES *FERRITIN	1220	3 mL (1.5 mL min.) serum from 1 SST AND 3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. Overnight fasting is preferred.
PV135	Immunoglobulin Profile	1000	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.
PV105	GLUCOSE FASTING (FBS) & GLUCOSE POSTPRANDIAL (PPBS)	180	2 mL (1 mL min.) plasma each for Fasting and Postprandial specimens from 1 Grey Top (Sodium Fluoride) tube. Overnight fasting mandatory. Collect fasting specimen. Note the time patient starts the meal. Collect PP sample exactly after 2 hours. Ship refrigerated or frozen
LIVE STOCK HERD HEALTH SCREENING TEST			
PV214	BOVINE BRUCELLOSIS	570	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.
PV215	PARATUBERCULOSIS (JOHNE'S DISEASE)	870	2 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.
PV216	BOVINE TUBERCULOSIS	2400	5 mL whole blood from 1 Green Top (Heparin) tube. Ship at room temperature.
PV222	FOOT AND MOUTH DISEASE (FMD)	690	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.
PV232	INFECTIOUS BOVINE RHINOTRACHEITIS (IBR)		2 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.
EXPORT TESTING SERVICES (PET TRAVEL TO ABROAD)			
PV220	CANINE TICK BORNE DISEASE PCR PANEL (Babesia, Ehrlichia, Hepatozoon)	2000	3 mL (2 mL min) Whole Blood in 1 Lavender Top (EDTA). Ship refrigerated
PV111	BLOOD SMEAR EXAMINATION FOR BLOOD PARASITES	300	Peripheral whole blood or 2 Peripheral blood smears. Give brief clinical history.
PV203	LEISHMANIA ANTIBODY, IgG	1600	2 mL (1 mL min.) serum from 1 SST Ship refrigerated or frozen.
PV091	KNOTT'S HEARTWORM TEST FOR MICROFILARIA	400	3 mL (1 mL min) Whole Blood in 1 Lavender Top (EDTA). Ship refrigerated
PV093	HEARTWORM ANTIGEN TEST	700	3 mL Whole Blood, Serum or Plasma. Ship refrigerated
PV205	RABIES SEROLOGICAL TESTING (FAVN)	16500	5 mL blood in SST or 3 mL Serum collected aseptically

TEST CODE	TEST NAME	TEST PRICE	TEST SPECIMEN
LAB ANIMAL HEALTH MONITORING AND CLINICAL TRIAL TESTING			
PV230	PATHVETS MOUSE HEALTH MONITORING PANEL A SEROLOGY: MHV, Sendai, LCMV, MPV 1-5, TMEV, PVM, EDIM, MAV 1, M. pulmonis PCR: C. rodentium, C. kutscheri, Salmonella spp, K. pneumoniae, Ps. aeruginosa, Staph aureus, P. multocida, Strep pneumoniae, Helicobacter spp	Details as per request, please contact 011-49058408	Opti-Spot, Feces, Oral Swab
PV131	HISTOPATHOLOGY, SLIDE PREPARATION AND INTERPRETATION		Submit specimen in 10% buffered formalin. Formalin should cover the specimen completely. Ship at room temperature. Brief clinical history in Histology Requisition Form is mandatory.
	ELECTROLYTES, URINE *SODIUM * POTASSIUM * CHLORIDE		10 mL (5 mL min) aliquot urine in a sterile screw capped container.

PROPER BLOOD COLLECTION METHOD USING VACUTAINER

TUBE/COMMON USES	COLOUR CODE	MIX BY INVERTING
Blood Cultures bottle *Aerobic # Anaerobic		8 to 10 times 
Citrate Tube Coagulation study		5 times
Vacutainer SST Gel Separator Tube (Clot activator tube, with gel barrier) Biochemistry, Hormone Assay, Serology		5 times
Serum Plain Tube (No Additive) Biochemistry, Progesterone, Phenobarbitone		5 times
Heparin Tube Bovine TB IFN- γ ,		8-10 times
EDTA Tube (Whole blood) CBC, Blood Parasite, ESR, PCR for Blood Parasites, Ammonia		8 to 10 times
Fluoride (glucose) Tube (Sodium fluoride is the antiglycolytic agent that will produce plasma) Used for Glucose testing		8 to 10 times

Route to Veterinary excellence



For all protocols involved in collection of Sample/Transportation and documentation, please refer to the contact number 011-49058408.



Ground Floor 40/50, Block B, Chittaranjan Park, Near Kotak Mahindra Bank,
New Delhi – 110019, Lab No. - 011-49058408,
E-mail: lpl.pathvets@lalpathlabs.com, Web: www.pathvets.in

TAT-It refers to the sample reaching PathVets Lab, CR Park, New Delhi.
Specimens, Method, Reporting time may change. Please contact customer care before sending the samples to avoid any inconvenience.